


Absorption Spectra of CuSO_4 and CoSO_4

Graph B



PART C. SPECTROPHOTOMETRIC ANALYSIS FOR COPPER

Sample Number #5 #6

Mass of paper plus compound _____

Mass of paper alone _____

Mass of compound taken _____

IMPORTANT! Concentration of standard Cu^{2+} solution _____ M

Table C1:


Sample No.	mL standard Cu^{2+}	mL 1 M HNO_3	Actual Cu^{2+} Concentration (M)
1	0.0	10.0	0.0
2	4.0	6.0	
3	7.0	3.0	
4	10.0	0.0	

Table C2: Observed Absorbance at λ_{max} (Use $\lambda = 645$ nm for the Cu^{2+} determination)

Sample No.	Trial 1	Trial 2	Trial 3	Average
1 (Dilute HNO_3)	0.0	0.0	0.0	0.0
2				
3				
4 (Standard Cu^{2+} solution)				
5 (Unknown)				
6 (Unknown)				

Prepare your calibration plot on the paper provided in the laboratory and draw the best straight line through the points for samples 1, 2, 3, and 4. Use the computer to prepare the same graph. Plot the data as illustrated here.

See page E-5.



Do your points fall on a straight line? If not, comment on the possible reasons for this? Do you need to redo some experimental measurements to correct the problem?

From the computer plot, record the following information:

• Slope of the line = $(\Delta \text{Absorbance} / \Delta \text{Concentration}) = \underline{\hspace{2cm}}$

• Intercept = $\underline{\hspace{2cm}}$ (This should be about 0)

(Note: This same information can also be obtained readily from the graph done on paper.)

The copper concentrations in your unknown samples (5 and 6) can be determined in two ways:

- You can read the concentrations from your graph.
- You can calculate them from the equation $y = mx + b$, where x is the unknown concentration, m is the slope of the line, y is the measured absorbance of the sample, and b is the intercept of the line.

Unknown Concentrations

Sample 5: $\text{Cu}^{2+} = \underline{\hspace{1cm}}$ M

Sample 6: $\text{Cu}^{2+} = \underline{\hspace{1cm}}$ M

How did you determine the results above? From your graph or from line slope? If from the line slope, show a sample calculation.

Calculating the Amount of Copper in the Unknown: Sample No. = $\underline{\hspace{2cm}}$

Grams of Cu^{2+} in your samples = (concentration)(0.01 liters of soln.)(63.54 g/mol)

Sample 5: grams of $\text{Cu}^{2+} \underline{\hspace{1cm}}$

Sample 6: grams of $\text{Cu}^{2+} \underline{\hspace{1cm}}$

Sample 5: wt. % $\text{Cu}^{2+} \underline{\hspace{1cm}}$

Sample 6: wt. % $\text{Cu}^{2+} \underline{\hspace{1cm}}$

Average % copper in the sample = $\underline{\hspace{2cm}}$

