

Name: _____ Section: _____ Date: _____

Mass of aspirin tablet= _____ g = _____ mg

TABLE 1

TUBE	CONTENTS	DILUTION FACTOR	CONCENTRATION (mg/ml)	ABSORBANCE AT 530nm
1	3mL FeCl ₃ solution	-	0.00	
2	1mL ASA standard 9mL FeCl ₃ solution	-	0.16	
3	5mL Tube #2 5mL FeCl ₃ solution	-	0.08	
4	5mL Tube #3 5mL FeCl ₃ solution	-	0.04	
5	5mL Tube #4 5mL FeCl ₃ solution	-	0.02	
A	1mL aspirin solution 9mL FeCl ₃ solution	10	?*	
B	0.5ml aspirin solution 9.5mL FeCl ₃ solution	20	?*	

*** NOTE**

The concentration of colored product in tubes A and B will be determined by comparing the absorbance with the absorbance of the "standard", or known samples (tubes 1-5), using a Standard Curve (see graph on following page). Tube B is twice as dilute as Tube A, so it should absorb half as much light.

DATA ANALYSIS- STANDARD CURVE

In order to determine the amount of ASA in the aspirin tablet, you will have to compare the absorbance of your aspirin samples to that of known concentrations of ASA (tubes 1-5). Plot the absorbance versus the concentration for tubes 1-5 on the grid below (concentration on x-axis, absorbance on the y-axis). Draw a straight “best fit” line through these points (a best fit line comes as close as possible, but *may not go through all data points*). This is your “standard curve”. It tells you the relationship between the absorbance of a salicylic acid solution and its concentration. Next, find the absorbance value for sample A on your graph. Based on your “standard curve”, read what the concentration value for that sample should be and record in Table 1.



